

Tuttnauer

LEADING IN
STERILIZATION
SINCE
1925

T-CLASSIC

9-Litre Portable Autoclave

Making your process easy, quick and affordable.

The T-Classic One, is a lightweight, portable steam autoclave designed as a cost-effective sterilization solution. Fully automatic with an easy one touch operation, delivering short cycles.

Ideally designed to meet portability needs of a beauty clinic or academic facility.

The T-Classic One can also be used in combination with another autoclave or as a back up.

Tuttnauer Performance with the simplicity and affordability your small practice needs



T-Classic One

Instruments sterilized and ready quickly



The T-Classic One offers competitively quick unwrapped cycle times, in as short as 18 minutes*. The fast point-of-use sterilization makes your turnaround a priority, increasing efficiency, while eliminating bottlenecks.

*Cycle times may vary, depending on load & weight.



The T-Classic One is a simple one step operated device. With one touch, the T-Classic One will begin an automatic cycle until completion. An electronic detector will turn off the power if water level is insufficient. The perfect autoclave for limited space.



The T-Classic One is the lightest tabletop autoclave available, weighing only 5.35 kg allowing for easy portability, while requiring minimal counter space.

T-CLASSIC



The **T-Classic One** offers multiple benefits for users. The Autoclave has 9 litre capacity, with a 126°C sterilizing temperature. The internal furniture consists of a basket as standard, however can also be supplied with cassette holders for instruments.

Technical Specifications

Maximum Width Overall	340mm
Maximum Height Overall	335mm
Maximum net weight	5.35 kg
Chamber Diameter/Height	210mm/230mm
Chamber capacity	9 Litres
Maximum Instrument length	228mm
Maximum Load weight	3kg
Sterilisation Temperature	126 °C
Operating Pressure	1.6 Bar
Minimum Sterilisation time	11 minutes
Normal Cycle Time	22 minutes
Voltage	230 Volts
Wattage	1500 Watts
Frequency	50 – 60 Hertz
Furniture	Basket



Top view, easy access



Basket for easy loading and unloading



Easy to use control panel

Model	Inner Chamber Dimensions ø x D (mm)	Volume (Litres)	Shipping Weight (kg)	Shipping Dimensions WxDxH (mm)	Single Phase Power Supply
T-Classic One	210 x 230	9	8.5	400 x 400 x 620	220-240V, 50-60Hz 1500W, 10A

Tuttנauer Europe B.V.
 Hoeksteen 11, 4815 PR
 PO Box 7191, 4800 GD Breda
 The Netherlands
 Tel: +31 765 423 510
 Fax: +31 765 423 540

International Sales and Marketing
 Website: tuttנauer.com
 Contact us at <https://tuttנauer.com/contact>



LTR330-0128-A