



Tuttnauer Steri-Center



Steri-Center

Adjustable Height Packing Table

Steri-Center is an ergonomic height adjustable packing table specifically designed for sterilization preparation packing of instrument baskets used in CSSDs and medical centers.

Steri-Center is fitted with two solid lifting pillars made of anodised aluminium with stainless steel feet. Maximum table load is 150kg per pillar. When overloaded, the lifting pillars are blocked from moving.

Steri-Center has a simple and smart height control system that allows the table height to be manually adjusted between 800mm and 1300mm. With the built-in memory function three different table height positions can be stored. This allows the table height to be conveniently positioned for several people using the same table.



The Steri-Center worktop surface is hygienic, made of water-resistant Hexpl plate, fitted with a scratch-resistant HPL layer, brushed stainless steel, or Corian® surface. All materials are corrosion-resistant and resilient to common cleaning and disinfection products. The table material is UV and heat resistant.

Steri-Center is available in the following dimensions (WxD):

1800 x 800 mm

1800 x 1000 mm

2000 x 800 mm

2000 x 1000 mm



HPL layer



Stainless Steel



Corian®



Accessories

A wide range of additional features can be added to the Steri-Center.

Electrical Sockets & UTP Connections
For electricity, Ethernet and phone lines



Inspection Lamp
With magnifier and adjustable light



LCD Mount
Flex arm mount for the attachment of an LCD screen

Computer Mount
For placing a PC under the table



Upper Structure

The table can be fitted with an upper structure for the attachment of trays and baskets used for holding small materials and supplies that are needed while working at the Steri-Center. The upper structure can also be fitted with fluorescent lighting and light shade (optional).



Drawer Unit

- Provides a large storage space for small equipment and other materials
- The drawer unit holds 3 drawers each of dimension: 600x400x100 mm (WxDxH)



Sterile Transport Trolley

A modular transport trolley made of aluminum, for transporting sterile items from CSSDs to operating rooms (OR).

- Sterile Transport Trolley is available as a single or double column trolley
- Available shelf configurations:
 - D-type ISO 1-2-3: 400 mm x 600 mm (W x D)
 - W-type ISO 1-2: 600 mm x 400 mm (W x D)
- Equipped with runners for up to 11 basket positions per column that are positioned so that 100 mm or 200 mm high baskets can be stacked in the trolley:
 - Holds 11 baskets with ISO 400 x 600 x 100 mm dimensions
 - Holds 6 baskets with ISO 400 x 600 x 200 mm dimensions
- The loading capacity is 150 kg per runner pair
- 4 Castors (2 castors with brakes): Wheel diameter 125 mm (Polypropylene with non-stripping thermoplastic rubber tire)
- Total height of 1900 mm (including castors)



The runners tilt to enable easy loading of baskets/trays and are locked once loaded to prevent baskets/trays from falling out.



Stainless Steel Transfer Baskets

Store wrapped instrument baskets, protecting them from damage. The standard ISO sized baskets fit into the trolley runners.

Stainless Steel Instrument Baskets

A wide range of sizes of stainless steel instrument baskets that meet strict European quality and design standards.



The Sterile Transport Trolley can be fitted with an isolating cover to protect it from contamination in transit.

Your Sterilization & Infection Control Partners



Company Profile

For over 85 years, Tuttnauer's sterilization and infection control products have been trusted by hospitals, universities, research institutes, clinics and laboratories throughout the world. Supplying a range of top-quality products to over 100 countries, Tuttnauer has earned global recognition as a leader in sterilization and infection control.

Global Partnerships

At Tuttnauer we feel that business means people dealing with people. We pride ourselves on our reputation for having long-lasting relationships with our customers based on commitment and trust, spanning over decades and distances.

Our Flexibility is Your Advantage

Beyond our unmatched range of products, we also manage complete turnkey solutions, including planning, design and installation of equipment, as well as consultation and feasibility studies, for projects of all sizes.

More from Tuttnauer:

Featuring Tuttnauer's range of cleaning, disinfection and sterilization solutions



Large Capacity Sterilizers
for Medical Use



Washer disinfectors for
hospitals and laboratories



85 and 160 liter pre & post
vacuum autoclaves

Tuttnauer Steri-Center

Tuttnauer USA Co.
25 Power Drive,
Hauppauge, NY 11788
Tel: +800 624 5836, +631 737 4850 Fax: +631 737 1034
E-mail: info@tuttnauerUSA.com

Tuttnauer Europe b.v.
Hoeksteen 11, 4815 PR
P.O.B. 7191, 4800 GD Breda
The Netherlands
Tel: +31 (0) 765 423 510 Fax: +31 (0) 765 423 540
E-mail: info@tuttnauer.nl

Distributed by: